

## SURGICAL SUCTION TUBE

tapered; Ø 2.5mm; 14.0cm

**SKU:** 085-018-25

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

|                   |                               |
|-------------------|-------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>          |
| <b>TYPE</b>       | <a href="#">Suction Tubes</a> |
| <b>KLS Martin</b> | 18-502-14                     |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                               |
|-------------------|-------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>          |
| <b>TYPE</b>       | <a href="#">Suction Tubes</a> |
| <b>KLS Martin</b> | 18-502-14                     |