

## SURGICAL SUCTION TUBE

tapered; Ø 3.0mm; 14.0cm

**SKU:** 085-018-30

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-503-14

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-503-14