

SYRINGE WITH GLASS BARREL

Luer-Lock connection; 50 ml

SKU: 087-000-50

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
KLS Martin	18-910-50
Medicon	09.78.37

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
KLS Martin	18-910-50
Medicon	09.78.37