

TAYLOR TISSUE FORCEPS

bayonet shaped; 1X2 teeth; dissector end

17.5cm (7"); width of jaw 1.6 +/-0.2

SKU: 020-317-17

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-383-17
AUTHORS	WAUGH
Medicon	07.27.18
SPECIFICATIONS	bajonett , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-383-17
AUTHORS	WAUGH
Medicon	07.27.18

SPECIFICATIONS

[bajonett](#), [toothed](#)