

## TC ADSON DRESSING FORCEPS

straight; serrated; jaw 1,5+/-0.2

12.0cm (4.8")


**SKU:** 021-070-12

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD151
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-250-12
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>Medicon</b>	06.20.62
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD151

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-250-12
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>Medicon</b>	06.20.62
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>