

TC ADSON DRESSING FORCEPS

straight; serrated; jaw 1,5+/-0.2

12.0cm (4.8")


SKU: 021-070-12

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD151
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-250-12
AUTHORS	ADSON
Medicon	06.20.62
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD151

MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-250-12
AUTHORS	ADSON
Medicon	06.20.62
SPECIFICATIONS	serrated , straight