

TC BABCOCK INTESTINAL GRASPING FORCEPS

straight; serrated

20.0cm (8")





SKU: 160-180-20

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tissue Seizing Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	30-190-20
AUTHORS	BABCOCK
Medicon	50.10.70
SPECIFICATIONS	serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tissue Seizing Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	30-190-20

AUTHORS

[BABCOCK](#)

Medicon

50.10.70

SPECIFICATIONS

[serrated](#)