

## TC BARRAQUER MICRO NEEDLE HOLDER

straight; serrated; with catch

18.0cm (7.2")


**SKU:** 098-034-18

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro-Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">BARRAQUER</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">with catch</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro-Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">BARRAQUER</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">with catch</a>