

TC BOZEMANN NEEDLE HOLDER

s-shaped; straight; serrated

20.0cm (8")



SKU: 051-185-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM100
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-659-20
AUTHORS	BOZEMANN
Medicon	10.37.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM100

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-659-20
AUTHORS	BOZEMANN
Medicon	10.37.70
SPECIFICATIONS	serrated , straight