


TC CRILE-WOOD NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



SKU: 051-060-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM017
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-18
AUTHORS	CRILE , WOOD
Medicon	10.20.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM017

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-18
AUTHORS	CRILE , WOOD
Medicon	10.20.68
SPECIFICATIONS	serrated , straight