


TC CRILE-WOOD NEEDLE HOLDER

straight; serrated

20.0cm (8")



SKU: 051-060-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM019
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-638-20
AUTHORS	CRILE , WOOD
Medicon	10.20.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM019

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-638-20
AUTHORS	CRILE , WOOD
Medicon	10.20.70
SPECIFICATIONS	serrated , straight