

TC CRILE-WOOD NEEDLE HOLDER

straight; serrated

23.0cm (9.2")

SKU: 051-060-23/LH

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CRILE , WOOD
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CRILE , WOOD
SPECIFICATIONS	serrated , straight