


TC CRILE-WOOD NEEDLE HOLDER

straight; serrated

26.0cm (10.4")



SKU: 051-060-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM030
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CRILE , WOOD
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM030
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CRILE , WOOD

SPECIFICATIONS

[serrated](#), [straight](#)