

TC CRILE-WOOD NEEDLE HOLDER LH

straight; serrated; for left hand use

18.0cm (7.2")



SKU: 051-060-18/LH

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CRILE , WOOD
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CRILE , WOOD
SPECIFICATIONS	serrated , straight