


TC CRILE-WOOD NEEDLE HOLDER

straight; serrated

15.0cm (6")

SKU: 051-060-15

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM016
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-638-15
AUTHORS	CRILE , WOOD
Medicon	10.20.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM016

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-638-15
AUTHORS	CRILE , WOOD
Medicon	10.20.65
SPECIFICATIONS	serrated , straight