

## TC CUSHING-TAYLOR DRESSING FORCEPS

bayonet shaped; serrated

18.0cm (7.2")

**SKU:** 021-117-18

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-256-18
<b>AUTHORS</b>	<a href="#">WANGENSTEEN</a>
<b>Medicon</b>	07.26.68
<b>SPECIFICATIONS</b>	<a href="#">bajonett</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-256-18

## AUTHORS

[WANGENSTEEN](#)

## Medicon

07.26.68

## SPECIFICATIONS

[bajonett](#), [serrated](#)