


TC DE BAKEY NEEDLE HOLDER

straight; serrated

23.0cm (9.2")



SKU: 051-114-23

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM036
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-23
AUTHORS	DE BAKEY
Medicon	10.20.73
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM036

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-23
AUTHORS	DE BAKEY
Medicon	10.20.73
SPECIFICATIONS	serrated , straight