


TC DE BAKEY NEEDLE HOLDER

straight; serrated

26.0cm (10.4")



SKU: 051-114-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM037
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-26
AUTHORS	DE BAKEY
Medicon	10.20.76
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM037

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-26
AUTHORS	DE BAKEY
Medicon	10.20.76
SPECIFICATIONS	serrated , straight