


TC DE BAKEY NEEDLE HOLDER

straight; serrated

26.0cm (10.4")



SKU: 051-114-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM037 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-646-26 |
| AUTHORS | DE BAKEY |
| Medicon | 10.20.76 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM037 |

| | |
|-----------------------|---|
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-646-26 |
| AUTHORS | DE BAKEY |
| Medicon | 10.20.76 |
| SPECIFICATIONS | serrated , straight |