


TC DE BAKEY NEEDLE HOLDER

straight; serrated

37.0cm (14.8")



SKU: 051-114-37

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	DE BAKEY
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	DE BAKEY
SPECIFICATIONS	serrated , straight