


TC DE BAKEY NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



SKU: 051-114-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM035
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-18
AUTHORS	DE BAKEY
Medicon	10.24.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM035

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-646-18
AUTHORS	DE BAKEY
Medicon	10.24.68
SPECIFICATIONS	serrated , straight