

TC DERF NEEDLE HOLDER

straight; serrated

12.0cm (4.8")



SKU: 051-010-12

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM008
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-632-12
AUTHORS	DERF
Medicon	10.12.62
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM008

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-632-12
AUTHORS	DERF
Medicon	10.12.62
SPECIFICATIONS	serrated , straight