

## TC FINOCHIETTO NEEDLE HOLDER

angled; serrated

27.0cm (10.8")




**SKU:** 051-211-27

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM092
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-653-27
<b>AUTHORS</b>	<a href="#">FINOCHIETTO</a>
<b>Medicon</b>	10.43.77
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM092

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-653-27
<b>AUTHORS</b>	<a href="#">FINOCHIETTO</a>
<b>Medicon</b>	10.43.77
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">serrated</a>