

## TC FLAT NOSE PLIER DOUBLE CONVERTIBLE

fig. 2; parallel grip

18.0cm (7.2")


**SKU:** 190-741-18

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Flat Noise Pliers</a>
<b>Aesculap</b>	LX182
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	22-511-02



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Flat Noise Pliers</a>
<b>Aesculap</b>	LX182
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	22-511-02