

TC FLAT NOSE PLIER DOUBLE CONVERTIBLE

fig. 3; parallel grip

18.0cm (7.2")




SKU: 190-742-18

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Flat Noise Pliers |
| Aesculap | LX183 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 22-511-03 |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Flat Noise Pliers |
| Aesculap | LX183 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 22-511-03 |