

## TC GERALD DRESSING FORCEPS

straight; serrated

18.0cm (7.2")

**SKU:** 021-102-18

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">CUSHING</a> , <a href="#">TAYLOR</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>AUTHORS</b>	<a href="#">CUSHING</a> , <a href="#">TAYLOR</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>