

TC GORNEY FACELIFT SCISSORS

straight; blunt/blunt; 1 blade toothed

20.0cm (8")




SKU: 011-332-20

Categories: [Open Surgery](#), [Scissors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Face-Lift Scissors
MATERIAL	Tungsten Carbide(TC)
KLS Martin	11-970-19
AUTHORS	GORNEY
SPECIFICATIONS	blunt/blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Face-Lift Scissors
MATERIAL	Tungsten Carbide(TC)
KLS Martin	11-970-19
AUTHORS	GORNEY

SPECIFICATIONS

[blunt/blunt](#), [straight](#)