

## TC HALSEY NEEDLE HOLDER

straight; serrated, for left hand use

13.0cm (5.2")



**SKU:** 051-040-13/LH

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Needle Holders</a>                      |
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>                |
| <b>AUTHORS</b>        | <a href="#">HALSEY</a>                              |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Needle Holders</a>                      |
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>                |
| <b>AUTHORS</b>        | <a href="#">HALSEY</a>                              |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |