


TC HEANEY NEEDLE HOLDER

curved; serrated

26.0cm (10.4")



SKU: 051-141-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-643-26 |
| AUTHORS | HEANEY |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-643-26 |
| AUTHORS | HEANEY |

SPECIFICATIONS

[curved](#), [serrated](#)