


TC HEANEY NEEDLE HOLDER

curved; serrated

26.0cm (10.4")



SKU: 051-141-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-643-26
AUTHORS	HEANEY
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-643-26
AUTHORS	HEANEY

SPECIFICATIONS

[curved](#), [serrated](#)