


TC HEANEY NEEDLE HOLDER

curved; serrated

20.0cm (8")



SKU: 051-141-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM081
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-643-20
AUTHORS	HEANEY
Medicon	10.31.71
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM081

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-643-20
AUTHORS	HEANEY
Medicon	10.31.71
SPECIFICATIONS	curved , serrated