

TC HEGAR-BAUMGARTNER NEEDLE HOLDER

straight; serrated

14.0cm (5.6")



SKU: 051-048-14

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM015
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-634-14
AUTHORS	BAUMGARTNER , HEGAR
Medicon	10.14.64
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM015

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-634-14
AUTHORS	BAUMGARTNER , HEGAR
Medicon	10.14.64
SPECIFICATIONS	serrated , straight