


## TC HEGAR-BAUMGARTNER NEEDLE HOLDER

straight; serrated

12.5cm (5")



**SKU:** 051-048-12

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-634-13
<b>AUTHORS</b>	<a href="#">BAUMGARTNER</a> , <a href="#">HEGAR</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-634-13
<b>AUTHORS</b>	<a href="#">BAUMGARTNER</a> , <a href="#">HEGAR</a>

## SPECIFICATIONS

[serrated](#), [straight](#)