

TC HEGAR-BAUMGARTNER NEEDLE HOLDER

straight; serrated

12.5cm (5")



SKU: 051-048-12

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-634-13
AUTHORS	BAUMGARTNER , HEGAR
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-634-13
AUTHORS	BAUMGARTNER , HEGAR

SPECIFICATIONS

[serrated](#), [straight](#)