


TC HEGAR-VASCULAR NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



SKU: 051-090-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HEGAR, VASCULAR
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HEGAR, VASCULAR
SPECIFICATIONS	serrated , straight