


TC HEGAR-VASCULAR NEEDLE HOLDER

straight; serrated

14.0cm (5.6")



SKU: 051-090-14

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HEGAR, VASCULAR
Medicon	10.25.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HEGAR, VASCULAR
Medicon	10.25.65

SPECIFICATIONS

[serrated](#), [straight](#)