

## TC IRIS DELICATE OPERATING SCISSORS

straight; sharp/sharp

11.5cm (4.6")




**SKU:** 011-244-11

**Categories:** [Open Surgery](#), [Scissors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Operating Scissors</a>
<b>Aesculap</b>	DO208
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	11-982-11
<b>AUTHORS</b>	<a href="#">IRIS</a>
<b>Medicon</b>	02.20.61
<b>SPECIFICATIONS</b>	<a href="#">sharp/sharp</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Operating Scissors</a>
<b>Aesculap</b>	DO208

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	11-982-11
<b>AUTHORS</b>	<a href="#">IRIS</a>
<b>Medicon</b>	02.20.61
<b>SPECIFICATIONS</b>	<a href="#">sharp/sharp, straight</a>