

TC IRIS DELICATE OPERATING SCISSORS

straight; sharp/sharp

11.5cm (4.6")




SKU: 011-244-11

Categories: [Open Surgery](#), [Scissors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Operating Scissors
Aesculap	DO208
MATERIAL	Tungsten Carbide(TC)
KLS Martin	11-982-11
AUTHORS	IRIS
Medicon	02.20.61
SPECIFICATIONS	sharp/sharp , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Operating Scissors
Aesculap	DO208

MATERIAL	Tungsten Carbide(TC)
KLS Martin	11-982-11
AUTHORS	IRIS
Medicon	02.20.61
SPECIFICATIONS	sharp/sharp, straight