

TC LANGENBECK NEEDLE HOLDER

straight; serrated; without catch

18.0cm (7.2")

SKU: 051-360-18

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM191
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-692-18
AUTHORS	LANGENBECK
Medicon	11.36.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM191

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-692-18
AUTHORS	LANGENBECK
Medicon	11.36.68
SPECIFICATIONS	serrated , straight