

TC LICHTENBERG NEEDLE HOLDER

straight; serrated; with catch

20.0cm (8")

SKU: 051-366-20

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM132
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-686-20
AUTHORS	LICHTENBERG
Medicon	11.26.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM132

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-686-20
AUTHORS	LICHTENBERG
Medicon	11.26.70
SPECIFICATIONS	serrated , straight