

TC LOCKLIN GUM SCISSORS

angled; shanks straight; serrated

16.0cm (6.6")


SKU: 011-891-16

Categories: [Open Surgery](#), [Scissors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Gum Scissors
MATERIAL	Tungsten Carbide(TC)
AUTHORS	DEAN
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Gum Scissors
MATERIAL	Tungsten Carbide(TC)
AUTHORS	DEAN
SPECIFICATIONS	angled