

TC LOCKLIN GUM SCISSORS

angled; shanks curved; serrated

16.0cm (6.6")




SKU: 011-890-16

Categories: [Open Surgery](#), [Scissors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Gum Scissors
MATERIAL	Tungsten Carbide(TC)
AUTHORS	LOCKLIN
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Gum Scissors
MATERIAL	Tungsten Carbide(TC)
AUTHORS	LOCKLIN
SPECIFICATIONS	angled