


TC MARTINI NEEDLE HOLDER

straight; smooth

20.0cm (8")



SKU: 051-120-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	MARTINI
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	MARTINI
SPECIFICATIONS	straight