

TC MASSON NEEDLE HOLDER

straight; serrated

27.0cm (10.8")

SKU: 051-150-27

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|-----------------------------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM090 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-648-27 |
| AUTHORS | MASSON |
| Medicon | 10.40.77 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM090 |

| | |
|-----------------------|-----------------------------------------------------|
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-648-27 |
| AUTHORS | MASSON |
| Medicon | 10.40.77 |
| SPECIFICATIONS | serrated , straight |