

TC MASSON NEEDLE HOLDER

straight; serrated

27.0cm (10.8")

SKU: 051-150-27

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM090
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-648-27
AUTHORS	MASSON
Medicon	10.40.77
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM090

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-648-27
AUTHORS	MASSON
Medicon	10.40.77
SPECIFICATIONS	serrated , straight