

TC MATHIEU-KOCHER NEEDLE HOLDER

straight; serrated; with catch

17.0cm (6.8")

SKU: 051-358-17

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM150
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-674-17
AUTHORS	KOCHER , MATHIEU
Medicon	11.18.67
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM150

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-674-17
AUTHORS	KOCHER, MATHIEU
Medicon	11.18.67
SPECIFICATIONS	serrated , straight