

TC MATHIEU-KOCHER NEEDLE HOLDER

straight; serrated; with catch

20.0cm (8")

SKU: 051-358-20

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM153
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-674-20
AUTHORS	KOCHER , MATHIEU
Medicon	11.18.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM153

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-674-20
AUTHORS	KOCHER, MATHIEU
Medicon	11.18.70
SPECIFICATIONS	serrated , straight