

TC MATHIEU NEEDLE HOLDER

straight; serrated; with catch

17.0cm (6.8")

SKU: 051-350-17

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM366
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-676-17
AUTHORS	MATHIEU
Medicon	11.20.67
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM366

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-676-17
AUTHORS	MATHIEU
Medicon	11.20.67
SPECIFICATIONS	serrated , straight