


## TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



**SKU:** 051-080-18

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Needle Holders</a>                      |
| <b>Aesculap</b>       | BM066   |
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>                |
| <b>KLS Martin</b>     | 20-642-18   |
| <b>AUTHORS</b>        | <a href="#">HEGAR</a> , <a href="#">MAYO</a>        |
| <b>Medicon</b>        | 10.18.68  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                 |                                |
|-----------------|--------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>           |
| <b>TYPE</b>     | <a href="#">Needle Holders</a> |
| <b>Aesculap</b> | BM066                          |

|                       |                                      |
|-----------------------|--------------------------------------|
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a> |
| <b>KLS Martin</b>     | 20-642-18                            |
| <b>AUTHORS</b>        | <a href="#">HEGAR, MAYO</a>          |
| <b>Medicon</b>        | 10.18.68                             |
| <b>SPECIFICATIONS</b> | <a href="#">serrated, straight</a>   |