


TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

20.0cm (8")




SKU: 051-080-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM067
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-20
AUTHORS	HEGAR , MAYO
Medicon	10.18.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM067

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-20
AUTHORS	HEGAR, MAYO
Medicon	10.18.70
SPECIFICATIONS	serrated , straight